

18-22 Kasım 2015

21.11.2015

12:36:54


**JESSICA
SCHILDE**

3.Ulusal
Klinik Mikrobiyoloji
Kongresi-2015



www.klin2015.org

18-22 Kasım 2015
TMMOB TSKMİD
Sakarya - Adapazarı



IRIDICA -

A NEW KIND OF CERTAINTY

Klimud 2015 / Belek - Antalya

Dr. Jessica Schilde

Molecular Application Specialist, EMEA

CHOOSE TRANSFORMATION

Abbott

Company Confidential ©2014 Abbott



IRIDICA – A new kind of certainty

Complete SOLUTION

- An end-to-end solution from specimen to answer
- Sample ID tracking throughout the process
- Automated, simple results reporting

Comprehensive COVERAGE

- Unprecedented Coverage (bacteria, fungi, and viruses)
- 1,000 pathogens

Short TURN AROUND TIME

< 6 HOUR Time to First Result

Direct SAMPLES, different SAMPLE TYPES

- Direct Sample Testing, culture not required.
- Whole blood; Sterile fluids and tissues; BAL, ETA; Plasma

Reliable IDENTIFICATION

- Identifies one or more pathogens in poly-microbial infections
- Not affected by antimicrobial pre-treatment



THE IRIDICA TECHNOLOGY

CHOOSE TRANSFORMATION

Abbott


Company Confidential ©2014 Abbott



PCR/ESI-MS Technology

- The IRIDICA technology combines two Nobel-prize winning technologies: PCR and ESI/MS


Step 1- Nucleic Acid Extraction from Direct Specimens: Blood, BAL, ETA, Tissues




Step 4- Signal Analysis by Triangulation and Database Matching to Identify Pathogens



Step 2- Unique Primer Design for Amplification of Broad Groups of Pathogens



Step 3- Mass Spectrometry to Weigh the DNA




CHOOSE TRANSFORMATION

Abbott



The IRIDICA Instruments



CHOOSE TRANSFORMATION

Abbott



IRIDICA Software

Workstation Home Screen

#	Test Order	Bead Beating Tube	Process Tube	Elution Tube	Reagent Cartridge	Assay Strip
1	33454 IRIDICA Fungal Assay IVD (BFNGS1)	ID: 116251 Lot: 707311	0000116251	0000116251	ID: 116251 Lot: 707311 Exp: 2024-10-06	
2	44522 IRIDICA BAC LRT Assay IVD (PNEU71)	ID: 116253 Lot: 707311	0000116252	0000116252	ID: 116252 Lot: 707311 Exp: 2024-10-06	
3	33454 IRIDICA BAC LRT Assay IVD (PNEU71)	ID: 693369 Lot: 897064	0000116253	0000116253	ID: 116253 Lot: 707311 Exp: 2024-10-06	
4						
5						
6						

Sample Prep Loading Screen

CHOOSE TRANSFORMATION

Abbott


Company Confidential ©2014 Abbott




The IRIDICA Instruments

Pre-Amplification

Bead Beater (BB)




Sample Prep (SP)



Post-Amplification

Thermal Cycler (TC)




Desalter (DS)

Private LAN



Analysis
Computer
(AC)



Mass Spec (MS)

CHOOSE TRANSFORMATION

Abbott




Pre-Amplification Workflow



Nucleic Acid Extraction and Assay Strip Setup

Pre-Amplification

Bead Beater (BB)



Sample Prep (SP)



CHOOSE TRANSFORMATION


Abbott

Company Confidential ©2014 Abbott




Post-Amplification Workflow

Amplification, Desalting and Analysis



Post-Amplification




CHOOSE TRANSFORMATION

Abbott



Organism Identification




CHOOSE TRANSFORMATION


Abbott



Broad Amplification

Primers are designed in a way that they span two universally conserved regions not too far apart from each other, but where the nucleotides in between the conserved regions are variable and differ with each unique species.







Abbott



Mass Spectrometry of PCR Amplicons, Signal Processing and Organism Identification by Triangulation



Triangulation Across Multiple Primers

Mass	Base Comp.	Quantity
30831.72	A27G30C21T21	40341
36378.33	A30G29C30T29	36030
-	-	No detection
33427.45	A29G30C25T24	40541
-	-	No detection
31271.05	A26G30C25T20	34377
24373.61	A16G23C21T19	8548
38952.31	A43G28C19T35	6305

Triangulation



S. aureus 6/6 primers
Quantity: 11558



SUMMARY REPORT



Sample ID:
test_4neg

Sample Type:	Negative Control	Assay Strip:	BCAL51201606179935D45541200001625
Entry Date:	23-Oct-2014 8:29 AM	DS Plate:	01351119900001192 (Position 2)
Spray Date:	23-Oct-2014 5:48 PM	Protocol Name:	IRIDICA BAC BSI Assay IVD
Result Date:	23-Oct-2014 5:57 PM	Protocol:	NPBCAL51 - 2.0.0.10349
Mass Spec ID:	MassSpec 1 (MS01)	Assay:	BCAL51 - V.48.50.42
Desalter ID:	Desalter 1 (DS01)	SPS Version:	SPS.2.0.10354
		GenX Version:	REL-GenX-V09R003-10213

IRIDICA Run Detail
information

Flags

WARNING: Negative Control sample has detections.

Assay Flags or Alerts

Bacteria

Detected Organism	Q Score	Level
Staphylococcus aureus	0.96	2

Bacteria and Fungi
Detection sections

Fungi

Detected Organism	Q Score	Level
Not Detected	—	—

Marker

Detected Marker	Result	Q Score	Level
KPC	n/a	—	—
mecA	Not Detected	—	—
vanA	n/a	—	—
vanB	n/a	—	—

Drug Resistance
Marker section

Control

Detection	Q Score	Level
Extraction Control	0.97	187

Extraction Control
Detection section

CHOOSE TRANSFORMATION

Abbott



Broad Assay Menu and Sample Types

IRIDICA ASSAY

BAC BSI

1

Blood stream infection

BAC SFT

2

Sterile fluids and tissues

BAC LRT

3

Lower respiratory tract

Fungal

4

Fungal

Viral IC

5

*Viral-
Immunocompromised*

Assay Type

Coverage

Selected Diseases/Patient groups

780 Bacteria and *Candida*,
4 Antibiotic Resistance Markers:
mecA, *vanA*, *vanB* and *kpc*

Identical coverage with
semi-quantitative threshold

> 200 fungi and yeast

13 viral reporting groups
covering > 130 viral species

Sepsis

Joint prosthetics infection

Pneumonia, IC patients

Pneumonia, IC Patients


IC Patients

CHOOSE TRANSFORMATION

Abbott



Parallel processing of different Assays & Samples




CHOOSE TRANSFORMATION

Abbott



IRIDICA Hands-on & Walk-away Times

**Bead Beater
(BB)**



**Sample Prep
(SP)**




**Thermal Cycler
(TC)**



**Desalter
(DS)**



**Mass Spectrometer
(MS)**



Hands on time (min):

15	10	1	5	1
-----------	-----------	----------	----------	----------

Walk away time (min):

5	110	140 – 200	35	30*
----------	------------	------------------	-----------	------------

*Hands on time and walk away time based on in-house data // * To first sample result*



STUDIES & PUBLICATIONS



The RADICAL Study

Rapid Diagnosis of Infection in the Critically Ill, a Multicenter Study of Molecular Detection in Bloodstream Infections, Pneumonia, and Sterile Site Infections*

Jean-Louis Vincent, MD, PhD, FCCM¹; David Brealey, MD²; Nicolas Libert, MD³; Nour Elhouda Abidi, MD⁴; Michael O'Dwyer, MD⁵; Kai Zacharowski, MD⁶; Malgorzata Mikaszewska-Sokolewicz, MD⁷; Jacques Schrenzel, MD⁸; François Simon, MD⁹; Mark Wilks, PhD⁵; Marcus Picard-Maureau, PhD¹⁰; Donald B. Chalfin, MD, MPH¹¹; David J. Ecker, PhD¹¹; Rangarajan Sampath, PhD¹¹; Mervyn Singer, MD²; the Rapid Diagnosis of Infections in the Critically Ill Team

Critical Care Medicine, 2015 Nov;43(11):2283-91


- A Multicenter Observational Trial to Evaluate the Potential Clinical Impact of IRIDICA
- Nine ICUs in six European countries as study sites
- Analysis of 616 bloodstream infection, 185 pneumonia and 110 sterile fluid and tissue specimens from 529 patients



The RADICAL Study Key Findings

Performance against culture in RADICAL:

- Sensitivity: >81%
- NPV: 97% for blood stream infections
- 3X greater recovery of pathogens





IRIDICA Early Assessment Program

M-H
Hannover Medical School

Comparison of the new PCR/ESI-MS platform IRIDICA with quantitative culture for detection of bacterial pathogens in bronchoalveolar lavage fluids of patients with suspected pneumonia

Philipp Kirschner¹, Sabrina Woltemate¹, Ines Yang¹, Stefan Ziesing¹, Tobias Weiler¹, Sebastian Sauerbaum¹

¹Institute of Medical Microbiology and Hospital Epidemiology and ²Department of Pulmonary Medicine, MHH Hannover Medical School, Carl-Neuberg-Str. 1, 30625 Hannover, Germany

ABSTRACT The Abbott PORE-MS (IRIDICA) assay (Fig. 1) is a qualitative in vitro test for the detection and identification of bacterial and *Candida* nucleic acids through PCR amplification and subsequent electrospray ionization mass spectrometry analysis. The test can also detect the presence of genes encoding resistance to certain antibiotics (medA, vanA, vanB) and *KPC* in association with specific bacterial detections.

The aim of this study was to evaluate IRIDICA for molecular diagnosis of respiratory pathogens. 115 bronchoalveolar lavage fluids (BAL) of patients with suspected pneumonia were tested. The results available after 1 hour were compared to the quantitative culture-based standard methods for pneumonia diagnosis. From a subset of 40 BAL, microbe composition was also performed by partial amplification of 16S rDNA and deep sequencing of 16S amplicons (HT sequencing) (Fig. 2).

Results

Of 115 patients, results were used to diagnose treatment of 105 (91%) for pneumonia patients (Fig. 1). The results were compared to the quantitative culture-based standard methods for pneumonia diagnosis. From a subset of 40 BAL, microbe composition was also performed by partial amplification of 16S rDNA and deep sequencing of 16S amplicons (HT sequencing) (Fig. 2).



RESEARCH ARTICLE

Evaluation of the Broad-Range PCR/ESI-MS Technology in Molecular D Infections

Elena Jordana-Lluch^{1,2,5}
Clara Marcó¹, M^a Jesús I
Vicente Ausina^{1,2,5}

CLINICAL BENEFITS OF RAPID PATHOGEN TESTING WITH PCR/ESI-MS

Dr. Mark Wicks Clinical Scientist, Microbiology at Barts Health NHS Trust in London, UK, talks about their experience of using PCR/ESI-MS technology over a period of 18 months. During this time for the IRIDICA study, the department also ran clinical samples of interest through the technology.

Which patient groups could potentially benefit from the PCR/ESI-MS technology?
There are a number of distinct clinical groups for which this technology promises to be quite rewarding, including patients with severe sepsis, pneumonia and compromised immune systems. Patients can be immunocompromised because they have recently had a transplant, or they could be haematology oncology or HIV patients. All the immunocompromised groups tend to be infected with unusual bacteria and fungi, which you might not normally look for. In addition, ordinary bacteria, which do not harm immunocompetent people, can have serious consequences in this group.

How would you summarise your experience with this technology?
In general it's been quite an exciting process and one which has caused a lot of interest in microbiology and in the different clinical departments in the hospital. In some cases it has been quite difficult to interpret the results, because there has been no technology like this before, so we have no framework with which to have our understanding of the results. It is a steep learning curve. Occasionally we have been baffled by an unexpected organism, one that is quite hard to culture in the laboratory usually.

How does this technology differ to conventional testing and what are its advantages?
There are a number of ways in which it differs from conventional microbiology testing. First is the speed of testing. We are used to a kind of 'guessing' approach, where nothing much happens for a minimum of 18 hours or perhaps 1 or 2, whereas with PCR/ESI-MS technology results are available in 6 hours. Another difference is that a lot of bacteria and especially fungi are very difficult to grow and are very slow growing. So with PCR/ESI-MS technology we

are getting a lot more positives coming through. With this technology you do not need to try and think of the names of an organism and try to grow it. You rely on the fact that this technology has a very broad coverage and therefore does the thinking for you. You just look for any infections at all.

How would you recommend using this technology to rule in or rule out infections?
At the moment it is too early to give clear guidelines. Obviously if you get the sample through

and you get a positive result then it's up to you to decide whether to act on it or not, or with any other test. This is a relatively easy comparison to rule out infections when you are relying on the high negative predictive value of the technology to rule out infections. This requires a lot of confidence and experience for people to act on the result and therefore to stop treating and maybe stop looking for further agents.

Why is the high negative predictive value so important in ruling out infections in patients?
The main hope is that we'll have enough confidence in the results to rule out infection. For example, a lot of patients in ICU, who were thought to be septic, actually don't have any infection at all. They might have SIRS, but that could be nothing to do with infection - it could be due to surgery, trauma or another reason.

But obviously the possibility of infection has to be considered and treatment may well be started. And then may be no underlying infection at all. What we hope is that our experience so far with ruling out infection will be maintained, and that we will have increasing confidence to act on the results and not to narrow antibiotic treatment or stop antibiotic treatment. In patient groups such as haematology oncology patients there is a huge amount of prophylactic anti-fungal treatment, despite the fact that clinicians don't really have any evidence of the patient having a fungal infection. However, the consequences of not treating an invasive fungal infection are so serious that they do not take the risk. This has implications for costs as well, if they can rule out having to treat these patients. Barts and London NHS Trust, for example, spends up to two million

pounds per year on antifungal treatment, much of which is almost certainly unnecessary. Another patient group where this high negative predictive value is important is preterm babies, who are admitted to neonatal intensive care units with possible sepsis. Often they are given five days or more antibiotic treatment. The consequences of unnecessary antibiotic treatment are extremely serious. It is not just the question of unnecessary treatments and encouraging antibiotic resistance. It can include the chances of getting necrotising enterocolitis and late-onset sepsis and death.

DISCLOSURE:
"Point-of-View" articles are part of the ICU Management Corporate Engagement Programme.

iridica. Identify early. Treat with confidence. Transform care.



IRIDICA - ENABLING EARLIER TRANSITION TO OPTIMIZED ANTIMICROBIAL THERAPY.

OVER 1,000 PATHOGENS, IN 6 HOURS, DIRECTLY FROM CLINICAL SAMPLES.

The mission of Abbott's IRIDICA is to create diagnostic solutions that can provide faster, more actionable results for critical infections. This is achieved by focusing on delivering an innovative approach to the detection and characterization of a broad array of microorganisms, contributing to Abbott's expanding role in molecular testing.

To learn more about IRIDICA, please visit <http://iridica.abbott.com>

CHOOSE TRANSFORMATION™



CHOOSE TRANSFORMATION

Abbott

Company Confidential ©2014 Abbott



Thank you!